

Material Used

- Natural high quality rubber latex powder Free

Indications:

- To cover hands during & protection in areas of potential risk associated with blood-borne pathogens and other substances or biohazards.

Contraindications:

- Not to be used in patients with known hypersensitivity to any of the materials used.

Instruction for Use:

- Wash your hands before donning gloves whenever possible. Pull a glove from the box with one hand. Most people start by holding the glove in their non- dominant hand.
- Holding the glove by the open edge with the fingers of your non-dominant hand, make sure that the glove is facing the right direction to facilitate the thumb without fumbling.
- Slide your opposite hand into the base of the glove, fingers spread slightly apart. Make sure that you find each opening with the tips of your fingers before you push your hand all the way in. Wriggle your fingers slightly to make sure that the glove is seated firmly and comfortably.
- Pull the other glove from the box and repeat the process with the other hand.
- Examine the gloves to make sure that there are no holes or tears prior to contact with a patient or any surface.

Care instructions:

- Single use only.

Conditions of Handling, Preservation and Storage:

- Not more than 5 cartoons on each other.
- Nice Ventilated place.
- Out of Sun light.
- Suitable temperature (18 – 32 oC).



Obelis^{SA}

53 Bd. Général Wahis
B-1030 Brussels, Belgium
Phone: 32.2.732.59.54
Fax: 32.2.732.60.03
E-mail: mail@obelis.net
www.obelis.net



See Instruction for Use



Non-Sterile

For Single Use Only

Contains Latex

Cautions

Manufacturer Address



Lot No.



Date of Manufacturing



Use by / Expiry date



Catalogue No.

ISSUE/REV DATE: 01 September 2016